

A SIMPLE IDEA TAKES SHAPE

Internal Bra Mesh Will Revolutionise Both Lift And Reduction Surgery

A South African invention is set to have significant global effect on how surgeons view the options available for breast surgery.

Designed specifically to support and lift the drooping breast without the need to use mammary prosthesis, the Breform mesh will also permanently counteract the undesirable recurrent Ptosis of the breast after Mastopexy or reductive surgery.

The idea of using a framework of sutures or mesh internally to stop ptosis is not a new one, but several attempts in the past have failed to deliver satisfactory results. The genius of the Breform lies not only in the choice of material used but also in taking the apparently simple step of pre-shaping the mesh for breast surgery ,and thus creating a genuine support structure.

Due to the shape of the Breform the results now achieved will not only lift but also expand the unwanted tubular appearance of the breast at the upper pole. Doctors using the Breform have noticed a very welcome additional advantage to the procedure that was not considered during the development stages, and this is the benefit of reduced scarring.

The support offered to the breast tissue by the Breform allows over time for the ongoing contraction, rather than broadening of procedural scars .

This action can also apply to the often noted problem of increase in Areola size. The reduced scarring and non expansion of the Areola result in aesthetically pleasing and youthful results for the long term.

Breform is a knitted multifilament Polyester mesh, of the kind used safely in surgery for over 40 years, it is also important to note that Breform will not interfere with diagnostic breast screening and is not palpable through breast tissue. Surgical techniques presently used for breast lift and reduction surgery need only minor alterations in order to safely and correctly position the Breform.


We do request that doctors wishing to utilise the Breform contact us in advance for information on the methodology of the procedure.



Before Breform Implant



After Breform Implant



TESTIMONIALS

"It is my belief that the Breform is the most significant development in lift surgery we have seen "

Dr Van Deventer - Cape Town (ISAPS presentation) S.Africa

"No recurrent ptosis observed, no serious complications encountered, the mesh cannot be felt in physical breast examination or seen on X-ray, excellent"

Dr H.P. de Bruijn - Enscheda Netherlands

"This system seems to be the solution for permanent mastopexy results that I have been looking for"

Dr S.Johannes - Johannesburg, S.Africa

"We have been waiting for a suitable internal mesh to use with Goes technique, with this internal bra system I think we have found it"

Prof.S. Monstray - Gent, Belgium

CLOVERLEAF MEDICAL 
SUPPLIERS OF SURGICAL PRODUCTS

WWW.CLOVERLEAFPRODUCTS.COM

TEL: 01494 876 990 FAX: 01494 876 991 EMAIL: MATT@CLOVERLEAFPRODUCTS.COM